



Human Growth Hormone Peptide CJC 1295 Without DAC For Body Performance Enhancement

Our Product Introduction

for more products please visit us on peptide-powder.com

Basic Information

- Place of Origin: China
- Brand Name: Hongbaiyi
- Certification: COA, HPLC MR
- Model Number: HBY-CJC 1295 Without DAC
- Minimum Order Quantity: 5 box
- Price: Negotiable
- Packaging Details: 5mg/vial, 10vials/box
- Delivery Time: 3-5 work days after your payment
- Payment Terms: , Western Union, T/T
- Supply Ability: 4000kg/Month



Product Specification

- Peptide Name: CJC 1295 Without DAC
- Molecular Weight: 3367.89
- Molecular Formula: C165H271N47O46
- CAS: 863288-34-0
- Specification: 5mg/vial, 10vials/box
- Appearance: White powder
- Storage: Cool Dry Place
- Shelf Life: 2 Years



More Images



Product Description

Human Growth Hormone Peptide CJC 1295 Without DAC For Body Performance Enhancement Basic Information Form Of CJC 1295 Without DAC

Product Name	CJC 1295 Without DAC
Cas No.	863288-34-0
Purity	99%min
Appearance	White powder
MF	C152H252N44O42
MW	3649.30
Shelf Life	2 years
Condition	Cool Dry Place

Brief Introduction Of CJC 1295 Without DAC

CJC-1295 is an injection peptide for increasing GH production. The peptide is a growth hormone releasing hormone (GHRH) mimetic or the like. In other words, it works the same way as GHRH, which can be called GHRH.

The primary use of CJC-1295 is to provide increased GH levels, which also leads to an increase in IGF-1 levels. These levels of increase can help reduce fat and, in some cases, also help muscle benefit. Typically, products in the GHRH category (including CJC-1295) are selected as alternatives to GH and are rarely used in conjunction with GH.

Another major GHRH product is Mod GRF 1-29, and in most cases I recommend CJC-1295. The duration of the product is different. Mod GRF has a roughly ideal short-acting time, allowing pulsed administration, whereas CJC-1295 has an extended duration of action that prevents this administration.

It is important to avoid confusing CJC-1295 with "CJC-1295 without DAC", which is not CJC-1295, but the wrong Mod GRF. When the peptide does not have a DAC, it is not CJC-1295.

CJC 1295 is sometimes referred to as "CJC-1295 and DAC", which is only CJC-1295.

Product Picture Of CJC 1295 Without DAC



When to use CJC-1295

The product is best suited for individuals who wish not to always inject and are seeking substantial support for GH production rather than maximum or near maximum increase. This is because it provides a flat blood level that does not match the pulsatile dose, which is required for maximum effect. However, a stable level can provide very good support for natural GH pulses.

In contrast, adverse side effects associated with GH overdose, such as nerve compression pain (such as carpal tunnel pain), excessive water retention, or reduced insulin sensitivity, compared to CJC-1295 use. The reason is that stimulating GH production is greater than the GH

production for individual situations. The solution is to stop using the problem until the problem is solved and reduce the dose when resuming use.

In order to continuously support GH production, the following recommended dose, CJC-1295 does not need to be recycled.

5. How to use CJC-1295

CJC1295 is usually set in vials containing 2 or 5 mg of lyophilized powder, although the amount may vary. The contents should be reconstituted by adding the appropriate amount of sterile or antibacterial water. If, for example, 2 mL is selected and the vial dose is 2 mg, the concentration of the resulting solution is 1 mg / mL or 1000 mcg / mL.

At the time of administration, the insulin syringe was used for stretching, and then the required amount was injected. In the above examples, a dose of 1000 mcg will require 1 mL of volume or & quot; 100 IU & quot; labeled on an insulin syringe.

Injection can be subcutaneous, intramuscular or intravenous according to personal preference. If necessary, if there is room, the peptide solution from other vials, such as vials of the GHRP product, may be inhaled into the same syringe. This reduces the total number of injections required.

When recommending CJC 1295, I usually recommend a dose of 1000 kg twice a week.

Package and Shipping of CJC 1295 Without DAC

Package: 5mg/vial, 10vials/box

Shipping: With 3-5 days after your payment



FAQ Of CJC 1295 Without DAC

Q1: Do you Accept Sample Order?

A: Yes. The smallest order we accept is one box.

Q2: Do you accept VISA business credit card?

A: Yes. VISA business credit card is in our payment.

Q3: Is there any discount?

A: Yes, if you want to purchase larger quantity, we can offer a better price for you.

Q4: How long does it take to me if I order?

A: About 3-5 days upon receipt of payment.

Q5: What kind of transportation do you have?

A: We have many ways of transportations, eg: Air, Sea, Special Lines, TNT, FEDEX, EMS, ect.



Shaanxi Hongbaiyi Biotech Co., Ltd.



18192109180



tracy@sxhongbaiyi.com



peptide-powder.com

Hengjia Business Building, No.115 Weiyang Road, E&T Development Zone, Xi'an, Shaanxi, China.